

# Lab 3

## Newton's Second Law

### Overview

In class, we related net force and change in momentum via the momentum principle:

$$\Delta \vec{p} = \vec{F}_{\text{net}} \Delta t \quad (1)$$

In the common case of constant mass  $m$ , this equation is the same as:

$$\vec{a} = \frac{\vec{F}_{\text{net}}}{m} \quad (2)$$

In today's lab, we will apply a constant force to cause a cart to accelerate down the track. We will measure acceleration, and investigate the relationship between the mass of the cart and the cart's acceleration.

### Hypothesize

#### Predictions

1. If we use the same force for each experiment but continually increase the cart's mass, what will we observe happen to the acceleration  $\vec{a}$ ?
2. If we use the same force for each experiment but continually increase the cart's mass, what will we observe happen to the change in momentum  $\Delta \vec{p}$  (assuming we measure momentum at the same time for each experiment)?

### Setup

1. Choose a hanging mass to act as the applied force throughout the experiment. A good force will move the cart across the track in about 2-3 seconds. The magnitude of your force is  $|\vec{F}| = mg$ , where  $g = 9.81 \text{ m s}^{-2}$ .
2. Record the mass of the cart (the mass of the cart by itself is 500 g. You will add additional mass throughout the experiment.)
3. Using the motion detector and Logger pro, measure the acceleration of the cart while the applied force pulls it down the track. Record the acceleration. Your instructor will show you how to use LoggerPro to measure acceleration.
4. Repeat the procedure with different cart masses. Try to obtain measurements for 5 different masses. Record your data in a table with one row per trial. The columns should be cart mass, measured acceleration, and net force.

### Analyze

Create a graph to visualize your data. On the x-axis should be the controlled variable (mass), on the y-axis should be the measured variable (acceleration). Use your data to calculate the constant force that was applied to the various cart masses. Compare your experimentally determined force to the theoretical force given by  $|\vec{F}| = mg$ .

### Questions

Revisit your predictions and compare them to your results.